

Intermittent Suction Unit Pediatric

DESCRIPTION AND APPLICATION

The Ohio Medical® Pediatric Intermittent Suction Unit (PED-ISU) is a medical suction device designed to deliver continuous or intermittent controlled vacuum from a hospital pipeline vacuum system. It provides the versatility of an intermittent mode for nasogastric drainage and a continuous mode for tracheal, pharyngeal, or surgical suctioning. The Pediatric model is designed to prevent vacuum adjustment above a preset vacuum limit level. This device is ideal for neonatal and pediatric patients or sensitive applications where limiting the peak vacuum level is important. The intermittent mode uses a pneumatic logic device to turn the suction “ON” and “OFF” automatically. This unique Unilogic timing module provides quiet and reliable controlled intermittent suction. The Pediatric model is adjustable from 0-150 mmHg, to meet most low clinical suction requirements. It also contains a positive pressure relief valve to reduce potential injury to the patient and equipment if momentarily attached to a positive pressure gas source.

The Pediatric Intermittent Suction Unit is a low vacuum/high flow regulator.

FEATURES AND BENEFITS

- **Baby Icon indicates Pediatric/Limited Model**
- **Vacuum Control**
 - Fully adjustable vacuum level from 0-150 mmHg
 - Peak vacuum limit adjustable from 80-150 mmHg, factory preset to 135 mmHg
- **Unilogic Timing Module**
 - Pneumatic, quiet, reliable operation
 - Returns to atmospheric pressure during “OFF” cycle
 - Timing cycles
 - Independently adjustable 3 to 30 seconds “ON” and “OFF” cycles
 - Factory preset at 15 seconds on/8 seconds off (begins in the off mode) Flow rates
- **Flow Rates**
 - (Intermittent) 0 to 13 L/min, preset at 8 L/min (ASTM)
 - (Continuous) 0 to 80 L/min (not adjustable)
- **Large, Easy to Read Gauge**
 - Color coded ranges
 - 5 mmHg increments
 - Glow-in-the-dark gauge increments and needle
 - Accuracy ± 5% full scale
- **Positive pressure and safety relief valve**
- **Modular Components**
 - Simplifies service procedures
 - Low cost replacement parts and kits
- **Durable Design and Construction**
 - Few mechanical parts
 - Rugged, shatter-resistant case
- **Ordering Options**
 - Available in various configurations to include adapters for connection to source vacuum. Also available with options for collection devices and overflow safety protection.



6700-1271-904
Pediatric Intermittent Suction Unit

TECHNICAL SPECIFICATIONS

Flow Rate

Intermittent – 0-13 L/min (preset at 8 L/min per ASTM)
 Continuous – 0-80 L/min (not adjustable) without fittings
 at full increase

Timing

3 to 30 seconds
 Independent on and off
 Preset 15 seconds on/8 seconds off
 Begins in the "OFF" mode

Vacuum Range

0-150 mmHg

Gauge Range

0-160 mmHg

Gauge Accuracy

±10 mmHg or ±5% Full Scale Deflection

Dimensions

6.6" H x 3.5" W x 4.8" D
 (16.8 x 8.9 x 12.2 cm)

Weight

1 lb 12oz (0.793 Kg)

F.O.B. Point

Gurnee, IL

Applied Standard

ISO 10079-3

* Specifications are nominal, subject to change without notice.

SHIPPING INFORMATION

Shipping Weight

2 lbs 8 oz (1.1 Kg) less fittings and adapters

Package Size

8" H x 3.5" W x 4.7" D
 (20.3 x 8.9 x 12 cm)


Packaging

One per carton

Terms

Net 30 days

Parts Configuration



© 2013 Ohio Medical Corporation.

This document contains information that is proprietary and confidential to Ohio Medical Corporation. Use of this information is under license from Ohio Medical Corporation. Any use other than that authorized by Ohio Medical Corporation is prohibited.

Ohio Medical, Ohio Medical Corporation and the Ohio Medical Corporation logo are registered trademarks of Ohio Medical Corporation.

Puritan Bennett is a registered trademark of Nellcor Puritan Bennett LLC.
 Chemetron is a registered trademark of Allied Healthcare Products, Inc.