


Product Specification



SPECIFICATIONS

Accurate Cracking Pressure
+/- 5 psi to 100 psi
+/- 10 psi from 100 psi to 200 psi
+/- 15 psi @ 200 psi and higher

Temperature Range
-100° F to +165° F
-73° C to +74° C

Pressure Settings
75, 125, 200, 250 and 375 psi

Orifice Size
5/16" Diameter

Flow Co-Efficient (CV Factor)
1.3 typ

Seal Material
Silicone

Part Numbers

232600	75 psi
232601	125 psi
232603	200 psi
232604	250 psi
232605	375 psi
232602	Pipe Away Adaptor
232606	Pipe Away Adaptor (1/4"F x 1/2"M use with 232601)

© 2012 Ohio Medical Corporation.

This document contains information that is proprietary and confidential to Ohio Medical Corporation. Use of this information is under license from Ohio Medical Corporation. Any use other than that authorized by Ohio Medical Corporation is prohibited. Ohio Medical Corporation, Ohio Medical and the Ohio Medical Corporation logo are registered trademarks of Ohio Medical Corporation.