

Reading the Digital Gauge

The **accuracy** of digital technology with the **familiarity** of a traditional analog gauge.

- Bars represent Low, Med, High and Full ranges
- Colors and bar graph are “quick glance” tools
- Easy to read from up to 30 feet away
- Large numbers to indicate **exact** suction level
- Color coding range on digital display corresponds to an analog gauge



NEONATAL (0-100mmHg)



LEVEL	RANGE	BAR DISPLAY
HIGH	85-100	8 BARS
HIGH	70-84	7 BARS
MED	55-69	6 BARS
MED	40-54	5 BARS
LOW	30-39	4 BARS
LOW	20-29	3 BARS
LOW	10-19	2 BARS
LOW	1-9	1 BAR
-----	0	NO BAR

PEDIATRIC (0-160mmHg)



LEVEL	RANGE	BAR DISPLAY
HIGH	140-160	8 BARS
HIGH	120-139	7 BARS
MED	100-119	6 BARS
MED	80-99	5 BARS
LOW	60-79	4 BARS
LOW	40-59	3 BARS
LOW	20-39	2 BARS
LOW	1-19	1 BAR
-----	0	NO BAR

ADULT (0-200mmHg)



LEVEL	RANGE	BAR DISPLAY
FULL	200	8 BARS & FULL
HIGH	160-199	7 BARS
HIGH	120-159	6 BARS
MED	100-119	5 BARS
MED	80-99	4 BARS
LOW	54-79	3 BARS
LOW	27-53	2 BARS
LOW	1-26	1 BAR
-----	0	NO BAR

SURGICAL (0-760mmHg)



LEVEL	RANGE	BAR DISPLAY
HIGH	700-760	8 BARS
HIGH	600-699	7 BARS
HIGH	500-599	6 BARS
HIGH	400-499	5 BARS
HIGH	300-399	4 BARS
HIGH	200-299	3 BARS
HIGH	100-199	2 BARS
HIGH	1-99	1 BAR
-----	0	NO BAR

© 2014 Ohio Medical Corporation. This document contains information that is proprietary and confidential to Ohio Medical Corporation. Use of this information is under license from Ohio Medical Corporation. Any use other than that authorized by Ohio Medical Corporation is prohibited. Ohio Medical Corporation and the Ohio Medical logo are registered trademarks of Ohio Medical Corporation. Amvex and the Amvex logo and Aeros are registered trademarks and Squire-Cogswell is a trademark of Ohio Medical Corporation. NFPA is a registered trademark of National Fire Protection Association Inc. ISO is a registered trademark of International Organization for Standardization.

TRUSTED BRANDS OF OHIO MEDICAL®



Ohio Medical Corporation
1111 Lakeside Drive - Gurnee, IL 60031-4099 USA
TollFree: 866-549-6446 - Fax: 847-855-6218
www.ohiomedical.com